



CLASSIC SOXHLET APPARATUS

Apparatus for classic, manual extraction using the Soxhlet method.

Tried and tested EFLAB quality for determining the fat content in food and feeds by means of solid/liquid extraction using an appropriate solvent.

Serial heating unit with 1-2-4-6 individually adjustable heaters. With interchangeable top moulds, airbath inserts, support rods and mains cables.

The cooling water supply system, Soxhlet glass parts and holders for fastening the glass parts to the support rods have to be ordered separately.

Opportunity to work with 250 and 500 ml glass flask
Optional Aluminum housing for round balloon use

SOXHLET GLASS SET, 4 AND 6 PIECES, 250 ML, CONSISTING OF:

- 1 Flat bottom flask 250 ml
- 1 Coiled-tube condenser
- 1 Extractor 150 ml

SOXHLET GLASS SET, 4 AND 6 PIECES, 500 ML, CONSISTING OF:

- 1 Flat bottom flask 500 ml
- 1 Coiled-tube condenser
- 1 Extractor 150 ml

